

The Self-contained Dental Treatment Unit P3 **Independence. Innovation. Class.**

The elegant P3 is a completely self-sufficient and mobile dental treatment unit. It requires no installation. Only 3 plugs need to be connected to the power supply for commissioning. All handpieces are freely selectable. Both the handpieces and the suction have full performance in continuous operation, like a fixed dental unit. Equipped with an integrated compressor and suction machine, the very powerful P3 offers comprehensive comfort due to its oil-free and extremely quiet operation (48 dB!). In addition, the P3's completely mobile carts allow maximum flexibility.

The Self-contained Dental Treatment Unit P3.
Independence. Innovation. Class.

