## Features

Integrated spray water system for syringe and instruments (check valve prevents the back flow of water)

1 Spray water tank: 0.75 Liter

Spray adjustment regulators for each instrument

Pressure regulator-filter with hose (6 × 4 mm) to air source

## Technical Data

Weight: 5.1 kg


Options: Electric scaler, curing light, electric micromotor, air motor, turbine, contra-angles, DENTA-TROLLEY, DENTA-CASE etc.

Denta-Boy 202

## The world's smallest dental unit

The DENTA-BOY 202 is quality in miniature.

Its small size means the DENTA-BOY 202 takes up very little space at your treatment location. It can be simply picked up by the handle and carried to the next treatment location. The DENTABOY 202 is equipped with a three-way syringe and two instruments of your choice.









Denta-Boy 202 Portable dental units