### Features

High quality stainless steel housing

1 Three-way syringe (air/water/spray) with sterilisable tips

Integrated spray water system for syringe and instruments (check valve prevents the back flow of water)

2 Spray water tanks: 2 × 0.75 litres

Instrument hoses detachable

Spray adjustment regulators for each instrument

Foot control with spray on/off function and chip

## Technical Data

Air consumption instruments: 35-55 I/min

Operating pressure handpieces: 2.5-3.5 bar

Each instrument is individually adjustable

Dimensions: 535mm × 580mm × 850mm

Options: Electric micromotor with light including display and 30 programs, electric scaler, peristaltic pump, turbine, air motor, curing light, tray etc.

## Denta-Cart 404

# Simply brilliant

The elegant DENTA-CART 404 treatment cart requires no installation and is totally mobile.

To use the DENTA-CART 404 simply connect it to an air source and plug it in the electrical socket! The DENTA-CART 404 simplifies your work and provides great comfort.

The instruments can be freely selected and are fully efficient even with continuous operation. The comfort is further enhanced by the great design and attention to detail.

Examples of this are the external removable instrument hoses and useful connections for external devices. Further the DENTA-CART 404 is equipped with two spray water tanks. If a spray water tank becomes empty, you can change to the second tank with a switch, continuing your treatment without interruption.



Mobile/Rollable

Denta-Cart 404

Mobile dental carts