Features

1 Three-way syringe (air/water/spray) with sterilisable tips

3 Instrument hoses with light (6-hole connector)

Integrated spray water system for syringe and instruments (check valve prevents the back flow of water)

1 Spray water tank: 0.75 litres

Spray adjustment regulators for each instrument

1 High-performance suction hose 300 l/min. Continuously adjustable. Unique BPR Swiss technology

Adapter for saliva ejector tips and oral-dry-cup

1 Suction fluid container: 1 litre (automatic overflow protection)

Pressure regulator-filter with hose (6 × 4 mm) for air connection

Sturdy, air- and water-tight plastic construction, will not scratch, dent or corrode

Detachable front cover

Easy-to-mount aluminium feet

Disc type foot control

Technical Data

Air consumption instruments: 35–55 I/min

Handpiece pressure: 2.5-3.5 bar

Each instrument is individually adjustable

Weight: 8 kg

Dilliensions: 370 jilli × 190 jilli × 500 jilli

Warranty: 2 years

Options: Electric scaler, curing light, electric micromotor, air motor, turbine, contra-angles, DENTA-TROLLEY, DENTA-CASE etc.

Denta-Port 303

The father of all portable dental units

The unmatched flexibility and quality of DENTA-PORT 303 are at the root of all BPR Swiss portable treatment units.

It can be equipped with all the instruments of your choice. Each instrument operates with full efficiency.

The DENTAPORT 303 is further equipped with the very powerful, worldwide unique BPR Swiss suction system (300 I/min!) and a three-way syringe. Activating the complete standard unit only requires plugging in the air connection.



Portable dental units Denta-Port 303