Features

The DENTA-DOCTOR STOOL 303 is set up or stowed away in less than 30 seconds

To comfortably transport, it can be stowed away in a convenient carrying bag

Height adjustment by means of a gas spring

Comfortable upholstery with premium covering

Upholstery available in many different colours

Technical Data

Weight: 5.65 kg

Diameter Seat: 340mm

Height adjustable: 540mm – 670mm

Warranty: 2 years

Options: Arm/Back rest

Denta-Doctor Stool 303

Simply sit perfectly everywhere

Sit comfortably. Disassemble. Transport.

The well conceived construction and high-quality workmanship combines robust construction, convenient transport, quick assembly and disassembly and high seating comfort with impressive lightness in a single dentist's chair.



Set up / disassembly in 30 seconds





Easy transport in corresponding carrying bag



Options: Back-/armrest

Contin

Continuous Operation



Plug-and-Play



Portable/Light

Denta-Doctor Stool 303