silent suction system and fresh water tanks

1 Three-way syringe (air/water/spray) with sterilizable tips

Unique BPR Swiss technology

1 High-performance suction hose 300 L/min with flexible cuspidor

2 Fresh water tanks: 2 × 2 litres

2 Suction fluid containers: 2 × 2 litres (automatic overflow protection)

Integrated suction with amalgam separator

Integrated oil-free compressor: Max. operating pressure of 4 bar

## Technical Data

Dimension: 355mm × 690mm × 895mm

Max. length of extended arm and cuspidor: 580mm

Warranty: 2 years

Options: Curing light

Master-Vac

## Worldwide unique assistant

The elegant MASTER-VAC is the first fully stand-alone suction cart whose functionality surpasses even fixed units.

The ceramic cuspidor, which is mounted on a movable arm, can be flexibly positioned under the patient's mouth. This makes rinsing the mouth easy for all patients, regardless of fitness, age or body size. During treatment, you as the dental professional will be assisted by a three-way syringe and an impressive suction capacity (300 I/min) that outperforms many fixed units. No costly installation work and still best comfort due to unique advantages.



Continuous Operation





Plug-and-Play



Mobile/Rollable

Mobile dental carts Master-Vac